RIDASCREEN[®] FAST Soya

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RIDASCREEN[®]FAST Soya

Brief information

RIDASCREEN[®]FAST Soya (Art. No. R7102) is a sandwich enzyme immunoassay for the quantitative analysis of soya proteins in untreated and processed food and beverages.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples)..... approx. 25 min
test implementation (incubation time)..... 30 min

Limit of detection: 0.31 mg/kg (ppm) soya protein

Limit of quantification: 2.5 mg/kg (ppm) soya protein

Specificity: The antibodies specifically detect soya proteins.
There is a cross-reactivity of 0.0017 % to legumes of the tribe Phaseoleae (various species of beans) and to legumes of the genus *Vicia faba* group (0.0003 %). To other legumes like peanut, lentil, pea, lupine as well as Milk proteins (casein, β -lactoglobulin) or egg proteins no cross-reactivity is measured.

1. Intended use

RIDASCREEN®FAST Soya (Art. No. R7102) is a sandwich enzyme immunoassay for the quantitative analysis of native and processed soya proteins in food. Such products are beverages as juice, wine, beer and food as sausages, dressings, bakery products, ice, chocolate, soup, sauce, margarine, etc.

2. General

Soya proteins can be present as an ingredient or as a contamination in raw and processed products. Several soya proteins can cause allergic reactions.

According to the EU Directive 2003/89/EG from November 10th 2003, soya and products thereof must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand.

Soybeans (mature seeds, raw) contain approx. 35 % proteins. Soya proteins have a similar amino acid composition as proteins from animal sources and have a high percentage of essential amino acids. Therefore, soya is one of the few plants which are considered as a valuable source of proteins. Thus, soya is widely used as substitute for proteins from animal sources (e.g. in tofu, soya milk or soya yoghurt). Besides the high protein content, soya beans are also rich in fat (approx. 20 %) and are used for the productions of oils and fats.

The usage of soya as food and feed has strongly increased in the last decades, leading to an increase in allergic reactions as well.

3. Test principle

The basis of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against soya proteins. By adding the standard or sample solution to the wells, present soya proteins will bind to the specific antibodies. The result is an antibody-antigen-complex.

In a washing step components not bound are removed. Then antibody conjugated to peroxidase (enzyme conjugate) is added. This antibody conjugate is bound to the antibody-antigen-complex. An antibody-antigen-antibody-complex (sandwich) is formed. Substrate/Chromogen solution is added after removing of any unbound enzyme conjugate in a washing step. Bound enzyme conjugate converts the chromogen into a blue product.

The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is proportional to the soya protein concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)
coated with anti-soya protein antibodies
- 5 x Soya protein standard solutions ^{*}), 1.3 ml each
0 mg/kg (ppm) (zero standard), 2.5 mg/kg, 5 mg/kg, 10 mg/kg, 20 mg/kg
soya protein in aqueous solution, ready to use
- 1 x Conjugate (0.7 ml)red cap
peroxidase conjugated antibody, 11x concentrate
- 1 x Red Chromogen Pro (10 ml) brown cap
Substrate/Chromogen solution, stained red
- 1 x Stop solution (14 ml)yellow cap
contains 1 N sulfuric acid
- 1 x Allergen Extraction buffer (125 ml)green cap
as a 10x concentrate
- 1 x Washing buffer (100 ml) brown cap
as a 10x concentrate
- 1 x Extraction Solution 2 (110 ml)..... blue cap
as a 2x concentrate

^{*}) The dilution factor of 100 for the sample has already been considered when labeling. Therefore, the soya protein concentration of samples can directly be read from the standard curve.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifugal vials
- shaker
- water bath (60 °C / 140 °F and 100 °C / 212 °F)
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- paper filter
- graduated pipette
- variable 20 - 200 µl and 200 - 1000 µl micropipettes
- if possible multichannel pipette and/or multistepper pipette

5.2. Reagents:

–distilled or deionized water

6. Warnings and precautions for the users

This test should be performed only by trained laboratory staff. The instructions for performing the test must be followed strictly.

The Extraction Solution 2 is harmful to health. It contains mercaptoethanol. It should be worked under a chemical hood and skin contact should be avoided (use gloves).

The stop solution contains 1 N sulfuric acid (R36/38, S2-26).

The washing buffer contains thimerosal. Avoid skin contact.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The red stained Substrate/Chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the Substrate/Chromogen solution prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 5

9. Preparation of Samples

The **Allergen Extraction buffer** is provided as a 10fold concentrate. Before dilution dissolve possibly formed crystals by heating (water bath 37 °C / 98.6 °F) and mix well. Dilute the heated concentrate 1:10 (1+9) with distilled water (e.g. 100 ml concentrate + 900 ml dist. water) or fill up the whole content of the bottle to 1250 ml. The diluted Allergen Extraction buffer can be used for approx. 12 weeks at 2 - 8 °C.

The **Extraction solution 2** is used **undiluted**.

9.1. Sample extraction

Unclean laboratory equipment can result in a contamination of samples or test components with soya. Therefore, gloves should be worn during the test procedure and the following recommendations should be applied:

- Thoroughly clean surfaces, glass vials, mixers and additional equipment before and after each sample
- Sample extraction and ELISA-procedure should be performed in separated rooms
- homogenize a representative amount of the sample (5 - 50 g)

Solid samples:

- weigh 1 g of sample, add 2.5 ml **undiluted** Extraction solution 2 and 17.5 ml preheated (60 °C / 140 °F) and finally diluted Allergen Extraction buffer, close the vial, mix vigorously and cook it for 10 min at 100 °C / 212 °F in a water bath
- let the sample cool down shortly (e.g. in an ice bath) and centrifuge 10 min at least 2500 g if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- dilute the extract 1:5 (1+4) with diluted Allergen Extraction buffer (e.g. 100 µl sample extract + 400 µl diluted Allergen Extraction buffer)
- pipette 100 µl extract per well

Liquid samples:

- take 1 ml of sample, add 2.5 ml **undiluted** Extraction solution 2 and 16.5 ml preheated (60 °C / 140 °F) and diluted Allergen Extraction buffer, close the vial, mix vigorously and cook it for 10 min at 100 °C / 212 °F in a water bath
- let the sample cool down shortly (e.g. in an ice bath) and centrifuge 10 min at least 2500 g if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)

- dilute the extract 1:5 (1+4) with diluted Allergen Extraction buffer (e.g. 100 µl sample extract + 400 µl diluted Allergen Extraction buffer)
- pipette 100 µl extract per well

The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 3 days. The extracts can be stored at -20 °C (-4 °F) for several months.

Remark:

Some sample matrices could influence the test result. Please do not hesitate to contact R-Biopharm AG if using e.g. buckwheat flour, spices or dark chocolate.

10. Test implementation

10.1. Test preparation for the ELISA

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **antibody enzyme conjugate** (bottle with red cap) is provided as a 11fold concentrate. Since the diluted enzyme conjugate solution has a limited stability, only the amount which actually is needed should be diluted. Before pipetting, the conjugate concentrate should be shaken carefully. For reconstitution, the conjugate concentrate is diluted 1:11 (1+10) with distilled water (e. g. 200 µl concentrate + 2 ml distilled water, sufficient for 2 microtiter strips).

The **washing buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml dist. water). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at a time.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard solution or prepared sample to separate duplicate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
4. Add 100 µl of the diluted enzyme conjugate to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl of the reddish Substrate/Chromogen solution to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win (Art. No. Z9999), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or soya contamination.

A further dilution and new detection of the samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 5. The samples should be diluted so that the results can be read in the linear part of the calibration curve.

Please note:

When working in accordance with the regulation stated, the dilution factor is 100. The soya protein concentration can be read directly from the standard curve - the sample dilution factor of 100 is already taken into account for the standard concentrations (see 4. *).

For sample dilutions of more than 1:100 the further dilution factor must be considered for the calculation of the soya protein concentration.

Remarks:

Samples tested negative still could contain a soya contamination below the limit of detection of the assay, or they might contain processed soya products like hydrolysates or proteinisolates for example which can contain protein fragments which are not further measurable in the sandwich format.

Recommendation:

In order to ensure a high analytical performance:

- each sample material should be analyzed in duplicates
- employ also soya protein free and soya protein containing (spiked) samples as test controls

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG