



Cultura[®] M

A small and handy incubator for dipslides

Cultura[®] M

Cultura M is a small, durable and reliable incubator suitable for the incubation of Orion Diagnostica's dip slides e.g. Uricult[®], Dentocult[®], Hygicult[®], Easicult[®].

Cultura M has been used for 10 years in the fields of medical and dental diagnostics, food bacteriology and hygiene monitoring, monitoring of coolants, lubricants and industrial liquids as well as industrial process.

Cultura M Incubator, Cat. No. 05866

Temperature range	25...45°C
Temperature stability	± 1°C
Voltage	230 V
Frequency	50 Hz
Power output	26 W / 0.2 A
Protective Class	II
Dimensions (mm)	external W: 310 H: 155 D: 168 internal W: 220 H: 120 D: 150
Case	Not easily flammable PC & ABS
Door	PC transparent
Thermometer	up to 60°C, in special tube
Base cover	ABS
Warranty	2 years from initial installation
Weight	1.1 kg



CE

This CE Mark indicates that this equipment complies with the following directive:

- EMC 89/336/EEC
- DIN EN 61010-4
- 61010-2010

